

## Small Chamber Sample Size and Shipping Instructions

One ASTM D6007 small chamber sample analysis requires three (3) specimens measuring 7-7/8" x 15" (200mm x 381mm).

### Shipping Instructions

- Use a marker to label each specimen with an identification code. Do not use stickers on test panels.
- Seal the edges of the samples with aluminum tape, if the samples arrive at the lab without edge sealing, they will be sealed with parafilm wax at the lab.
- All specimens must be bagged in a 6-mil plastic zip-lock bag and then appropriately boxed for shipping.
- All shipments must be pre-paid; no COD's.
- **Please notify the ITC before shipping any test material.**

**Please contact us at:**

labdir@cpamail.org  
(703) 724-1128 ext. 295

**Please send samples to:**

ITC  
73 Lawson Road, #101  
Leesburg, VA 20175 USA  
Attn: Lab Director